STORZ

KARL STORZ—ENDOSKOPE

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INSTRUCTION MANUAL
C-MAC® S IMAGER 8402 XS/8403 XS

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Thank you for your expression of confidence in the KARL STORZ name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern and high-quality instrument from KARL STORZ.

This manual is designed to serve as an aid in the proper handling, cleaning and sterilization of the C-MAC® S IMAGER. All required details and all actions on your part are clearly presented and explained.

We thus ask that you read this manual carefully before proceeding to work with the instrument. Keep this manual available for ready reference.

It is recommended to check the suitability of the product for the intended procedure prior to use.

KARL STORZ is constantly working on the further development of all products. Please appreciate that changes to the scope of supply in form, equipment and technology are possible for this reason. Therefore, no claims may be derived from the information, figures, and descriptions in this manual.

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1	Controls	4
2	Symbol description.....	6
2.1	Pictograms	6
2.2	Symbols for safety instructions	7
3	Normal use.....	8
3.1	Intended use.....	8
3.2	Contraindications.....	8
3.3	Target groups	8
3.4	User qualification	8
3.5	Training in the operation and function of the device	9
3.6	Patient population.....	9
3.7	User profile of physician and assistants	9
4	Safety Instructions	10
5	Handling.....	12
5.1	Possible combinations.....	12
5.2	Inserting the C-MAC® S IMAGER into the blade	12
5.3	Connecting the C-MAC® S IMAGER to the monitor.....	13
5.4	Connecting the power supply unit	14
5.5	Switching on the monitor.....	14
5.6	Disassembly of the C-MAC® S IMAGER	14
6	Reprocessing	15
6.1	Overview table	15
6.2	General warnings.....	16
6.3	Preparation for cleaning and disinfection	18
6.4	Manual wipe-down disinfection.....	18
6.5	Manual precleaning.....	18
6.5.1	No pretreatment with ultrasound	18
6.5.2	Manual cleaning.....	18
6.6	Manual disinfection	19
6.6.1	Low Level Disinfection (LLD)	19
6.6.2	Manual High Level Disinfection	19
6.7	Machine cleaning/ disinfection	20
6.8	Assembly, inspection and care	20
6.9	Packaging systems.....	21

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Contents

6.10	Sterilization	21
6.10.1	Chemical low-temperature sterilization with peracetic acid – STERIS® SYSTEM 1®	21
6.10.2	Chemical low-temperature sterilization with peracetic acid – STERIS® SYSTEM 1E®	22
6.11	Reprocessing limits.....	22
7	Service and repair	23
7.1	Servicing and repair.....	23
7.2	Limitation of liability.....	23
7.3	Warranty	23
7.4	Maintenance.....	23
7.5	Directive compliance.....	24
7.6	Standard compliance.....	24
7.7	Disposal.....	25
7.8	Repair program.....	25
7.9	Important information.....	25
8	Technical description	26
9	Appendix	27
9.1	Troubleshooting	27
9.2	Accessories	28
9.3	Information on electromagnetic compatibility (EMC)	29
10	Subsidiaries.....	42

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Controls

- ① Telescope
- ② C-MAC® S IMAGER for use with the video laryngoscope blades for single use
- ③ Extraction mechanism
- ④ Connecting cord for video monitor
- ⑤ Mount for the C-MAC® S IMAGER
- ⑥ C-MAC® S video laryngoscope blade for single use (not included in the scope of supply of the C-MAC® S IMAGER)

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Consult instructions for use!

Applied part type BF

Control of pollution caused by electronic information products (China RoHS)

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2.2 Symbols for safety instructions

WARNING: A Warning indicates that the personal safety of the patient or personnel may be endangered. Failure to observe a Warning could result in injury to the patient or personnel.

CAUTION: A Caution indicates that particular procedures or precautions must be followed to avoid possible damage to the instrument set.

NOTE: A Note indicates special information regarding operation or clarifies important information.

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3 Normal use

3.1 Intended use

The C-MAC® S IMAGER is used in conjunction with a C-MAC® S video laryngoscope blade for single use 051113-10/051114-10/051116-10 and the C-MAC® monitor 8403 ZX for endotracheal intubation and for inspection of the oropharynx.

3.2 Contraindications

No contraindications directly related to the product are currently known. The use of the C-MAC® S IMAGER in conjunction with a C-MAC® S video laryngoscope blade for single use 051113-10/051114-10/051116-10 and the C-MAC® monitor 8403 ZX or other compatible display devices is contraindicated if, in the opinion of the responsible physician, the patient's health is endangered through its use, e.g., due to the patient's general condition or if laryngoscopy as such is contraindicated.

3.3 Target groups

This instruction manual is intended for physicians, medical support staff and sterilization assistants.

3.4 User qualification

The C-MAC® S IMAGER and the accessories for single use operated in combination may only be used by persons with an appropriate medical qualification and who are acquainted with laryngoscopic techniques. The information given in this manual only serves to instruct in the correct handling, cleaning and sterilization of the C-MAC® S IMAGER. It is not intended as an introduction to the techniques of laryngoscopy. All users must have been instructed in the operation and use of the C-MAC® S IMAGER 8402 XS/8403 XS with the corresponding C-MAC® S video laryngoscope blade for single use, otherwise this medical device may not be used.

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The C-MAC® S IMAGER 8402 XS/8403 XS with the corresponding C-MAC® S video laryngoscope blade for single use must be used in accordance with the airway algorithm in the respective country of application or the official specifications for airway management.

3.5 Training in the operation and function of the device

Your local representative or responsible KARL STORZ member of staff is available to provide training and information on further training alternatives.

3.6 Patient population

Men and women as well as adolescents.
State of health: suitable for laryngoscopic interventions, in the physician's judgement.

3.7 User profile of physician and assistants

- Recognized medical skills in the field of application (specialist physician, qualified medical staff).
- Adequate powers of comprehension to rationally assess the surgical situation in hand.
- Adequate language skills in one of the languages used on the device and in the instruction manual.
- Thorough training in the operation of the instruments.
- Knowledge of the contents of the instruction manual.
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic).

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4 Safety Instructions

Risks of infection

- **Warning:** These instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties. Inspect instruments for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the instruments before initial use and before and after every subsequent use using validated procedures.
- **Warning:** Incorrectly reprocessed instruments expose patients, users and third parties to a risk of infection as well as the danger of an instrument malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.
- **Warning:** (if using sterile disposable products:) Before use, check the expiration date and ensure that the packaging is undamaged. If the expiration date has passed or the packaging is damaged, the product must not be used.
- **Warning:** (if using sterile disposable products:) The use of non-sterile products poses a risk of infection for patients. Do not remove the products from the sterile packaging until immediately before use.

Risks of injury

- **Warning:** Optical radiation poses a risk of injury to eyes. Never look directly into the light outlet at the blade tip or into the light outlet of the C-MAC® S IMAGER.
- **Warning:** The C-MAC® S IMAGER must not be used for medical applications without a C-MAC® S video laryngoscope blade for single use.
- **Warning:** The C-MAC® S IMAGER must be checked prior to each use to ensure that it functions correctly and has been properly cleaned, disinfected and/or sterilized. Defective instruments must be withdrawn from use and inadequately reprocessed instruments must be reprocessed again.

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Risks due to combination of system components

- Warning: Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product. All relevant instruction manuals must be read carefully and the instructions given must always be followed precisely. Check the functioning of the products used in combination.
- Warning: The instruction manuals and interface specifications for medical devices and/or system components used in combination must be observed precisely. Suitable system components are listed in Section 5.1 'Possible combinations'.
- Warning: Before using the instrument on the patient, a test must be carried out to ensure the proper functioning of the complete C-MAC® S IMAGER system, inserted in a C-MAC® S video laryngoscope blade for single use and connected to a C-MAC® monitor.

Damage to the instrument

- Caution: The C-MAC® S IMAGER must not be cleaned in an ultrasound bath.
- Caution: The C-MAC® S IMAGER must not be steam-sterilized.
- Caution: Do not place any objects on the product or the packaging containing the product.
- Caution: Do not store in the proximity of chemicals, disinfectants or radioactive emissions.

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* Software-Stand des Monitors 8403 ZX:

704v300 oder höher

* Software version of the monitor 8403 ZX:
704v300 or higher

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5.2 Inserting the C-MAC® S IMAGER into the blade

Open the packaging of the desired C-MAC® S video laryngoscope blade.

NOTE: The instruction manual of the C-MAC® S video laryngoscope blade for single use must be observed (Art. no. 96306060DF)!

Push the C-MAC® S IMAGER ② into the mount on the C-MAC® S video laryngoscope blade, until the IMAGER locks into place and you hear an audible 'click'.

NOTE: Check that the IMAGER is securely in place in the C-MAC® S video laryngoscope blade by gently moving it back and forth (to the right and left). The distal end of the C-MAC® S IMAGER must lie flat on the lens of the C-MAC® S video laryngoscope blade.

NOTE: The C-MAC® S IMAGER can be removed from the C-MAC® S video laryngoscope blade during operation and inserted into another C-MAC® S video laryngoscope blade.

WARNING: Danger of glare! Optical radiation poses a risk of injury to eyes. Never look directly into the light outlet at the blade tip or into the light outlet of the C-MAC® S IMAGER.

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i **NOTE:** The high-intensity light emitted by the C-MAC® S IMAGER can result in increased temperatures at the light outlet. However, these always stay below the value which could cause burns.

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5.3 Connecting the C-MAC® S IMAGER to the monitor

! **NOTE:** The instruction manual of the C-MAC® monitor must be observed (Art. no. 96076008D)!

To switch on the device, insert the connecting cord ④ of the C-MAC® S IMAGER into the corresponding connection socket ⑦ on the C-MAC® monitor 8403 ZX.

The connecting cord ④ can be connected or disconnected during operation.

! **CAUTION:** The C-MAC® S IMAGER 8403 XS can only be connected to the C-MAC® monitor 8403 ZX.

i **NOTE:** The C-MAC® S IMAGER 8402 XS can be connected to both the monitor 8402 ZX and the monitor 8403 ZX.

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The monitor 8403 ZX enables two CMOS video units to be connected with Lemo connectors.

i **NOTE:** All CMOS videoscopes from KARL STORZ can be connected to both connections. Video units from other manufacturers cannot be connected.

i **NOTE:** If two cameras are connected, the user must select the camera for which the images are to be displayed in the menu (refer to the instruction manual for the C-MAC® monitor 8403 ZX, Art. no. 96076008D).

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5.4 Connecting the power supply unit

Connect the cord of the power supply up to the connection socket on the C-MAC® monitor and plug it into the power supply. The power cord can be fixed to the connection socket using the union nut.

The battery is charged automatically as soon as the C-MAC® monitor is connected to the power supply. This is displayed by the lighting up of the orange-colored LED.

NOTE: The C-MAC® S IMAGER can also be used during the charging phase.

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5.5 Switching on the monitor

Switch on the monitor by pressing the power switch (the green pilot lamp lights up).

(For more information and further connection and control options of the monitor, refer to the instruction manual 'C-MAC® Monitor 8403 ZX', Art. no. 96076008D).

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5.6 Disassembly of the C-MAC® S IMAGER

The C-MAC® S IMAGER can be pulled out of the C-MAC® S video laryngoscope blade for single use by hooking a finger into the ring of the extraction mechanism ③.

! WARNING: Danger of glare! Optical radiation poses a risk of injury to eyes. Never look directly into the light outlet at the blade tip or into the light outlet of the C-MAC® S IMAGER.

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- This reprocessing step is part of the validated process
- This reprocessing step may be carried out additionally in order to aid cleaning
- X On the grounds of potential damage, this reprocessing step must NEVER be performed

No symbol: This reprocessing step has neither been validated in respect of efficacy nor in respect of material compatibility by KARL STORZ.

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6.2 General warnings

WARNING: Risk of infection: Incorrectly reprocessed medical devices expose patients, users and third parties to a risk of infection as well as the risk that the medical device may malfunction. Observe the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' and the accompanying documentation.

WARNING: Risk of infection: These medical devices are not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties. Inspect medical devices for visible contamination. Visible contamination is an indication that reprocessing has not been carried out or has been carried out incorrectly. Reprocess the medical devices before initial use and before and after every subsequent use using validated procedures.

WARNING: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety must be observed.

CAUTION: The C-MAC® S IMAGER has gold contacts which prevent corrosion. Particularly with manual reprocessing, however, it must be ensured that the contacts are rinsed thoroughly with microbiologically pure/sterile water during the final rinse. Following this, the contacts and the instrument must be dried carefully, preferably with sterile compressed air.

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Reprocessing

CAUTION: When preparing and using the solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life. Prolonged immersion and incorrect concentration may result in damage. Bear in mind the microbiological range of action of the chemicals used.

CAUTION: Danger of damage to medical devices: The use of chemicals which have not been approved by KARL STORZ may cause damage to the medical devices. Only use chemicals approved by KARL STORZ for reprocessing. A complete list can be found on the internet at www.karlstorz.com.

CAUTION: National laws and regulations must be observed.

NOTE: The manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments' can be downloaded or requested by visiting www.karlstorz.com.

NOTE: The connector of the C-MAC® S IMAGER 8402 XS/8403 XS is, in principle, completely waterproof. The connector is also resistant to the process chemicals used in the validated reprocessing methods employed. During the validated reprocessing methods, the connector can, in principle, remain unprotected.

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6.3 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

6.4 Manual wipe-down disinfection

Wipe-clean the exterior surfaces of the medical device with a disposable cloth moistened with disinfectant or a ready-to-use soaked disinfectant cloth. Due to their protein-fixating effect and possible material incompatibility, alcohol-based agents must not be used. The chemical manufacturer's specifications regarding material compatibility and necessary exposure time for the disinfectant must be observed. At the end of the necessary exposure time, wipe the surface with a dry low lint cloth.

6.5 Manual precleaning

6.5.1 No pretreatment with ultrasound

For technical reasons, the medical device is not suitable for ultrasound treatment.

6.5.2 Manual cleaning

The medical device must be completely immersed in a cleaning solution. Through disassembling and opening the instrument, and filling in a targeted manner, it must be ensured that even surfaces with restricted access and any lumina are covered and no air bubbles are present. After the necessary exposure time as specified by the manufacturer, clean the instrument mechanically with the aid of brushes or a sponge. Finally, it must be rinsed with cold water to ensure neutralization.

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6.6 Manual disinfection

The medical device must be completely immersed in a disinfectant solution. Through disassembling and opening the instrument, and filling in a targeted manner, it must be ensured that even surfaces with restricted access and any lumina are covered and no air bubbles are present. At the end of the necessary exposure time as specified by the manufacturer, the medical device must be rinsed several times in order to remove all chemical residues. To this end, water of the best possible quality must be used, taking into account the national regulations. Finally, all of the surfaces, joints, openings, channels and lumina are dried completely with (preferably sterile) compressed air in accordance with the national regulations. The cleaning gun with accessories (Art. no. 27660) is ideal for this purpose.

6.6.1 Low Level Disinfection (LLD)

KARL STORZ recommends using medical grade cleaning wipes with 0.1 % hypochlorite for low level disinfection of the C-MAC® S IMAGER. An example of such wipes are Clorox Healthcare Germicidal Wipes.

As a general rule, we recommend high level disinfection for the C-MAC® S IMAGER. However, if the C-MAC® S IMAGER is used, as intended, with a sterile C-MAC® blade for single use and no visible residues can be seen on the IMAGER after use, low level disinfection is acceptable.

6.6.2 Manual High Level Disinfection

- Cidex
- Cidex OPA
- Resert XL HLD

STERIS® SYSTEM 1E®

The following trays are compatible:

- C1220E Directed Flow Processing Container/Tray
- C1140E Flexible Endoscope Processing Container/Tray
- C1160E Universal Flexible Endoscope Tray

NOTE: The C-MAC® S IMAGER 8402 XS/8403 XS is not compatible with the C1200E General Processing Container/Tray.

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- Contenedor/bandeja C1220E Directed Flow Processing
- Contenedor/bandeja C1140E Flexible Endoscope Processing
- Bandeja C1160E Universal Flexible Endoscope

NOTA: El C-MAC® S IMAGER 8402 XS/8403 XS no es compatible con el contenedor/bandeja C1200E General Processing.

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Reprocessing

i NOTE: The Quick Connect Kit is not required as the instrument does not have any channels.

6.7 Machine cleaning/ disinfection

! **CAUTION:** Damage to the instrument: The C-MAC® S IMAGER may be cleaned and chemically disinfected in the washer and disinfectant at a maximum of 65 °C.

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

Machine cleaning/chemico-thermal disinfection

A maximum temperature of 65 °C must not be exceeded. The relevant national requirements and the chemical manufacturer's instructions must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

i NOTE: If necessary, the instrument must be dried off afterwards by hand.

6.8 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Clean the optical surfaces of the image chip and LED illuminator on the C-MAC® S IMAGER with a cotton tip applicator soaked in 70 % isopropyl alcohol.

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Reprocessing

- Damaged or corroded medical devices must be withdrawn from use.
- Afterwards, a functional check must be carried out.

i NOTE: The electrical contacts in particular must be thoroughly dried.

i NOTE: Use items from the catalog 'Care, Sterilization and Storage Techniques' during care procedures.

6.9 Packaging systems

Only standardized and approved packaging materials or systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

6.10 Sterilization

! **CAUTION:** Damage to the instrument: The C-MAC® S IMAGER 8402 XS/8403 XS must not be steam-sterilized.

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The procedure must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

6.10.1 Chemical low-temperature sterilization with peracetic acid – STERIS® SYSTEM 1®

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i NOTE: Sterilization is not possible on surfaces which have been greased and lubricated.

i NOTE: The C-MAC® S IMAGER 8402 XS/8403 XS is not compatible with the C1200E General Processing Container/Tray.

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6.10.2 Chemical low-temperature sterilization with peracetic acid – STERIS® SYSTEM 1E®

More detailed information on the selection of the sterilization parameters and the STERIS® Quick Connect Kit (QCK) required to rinse out the lumina is available from the manufacturer STERIS®.

-  **NOTE:** Sterilization is not possible on surfaces which have been greased and lubricated.
-  **NOTE:** The C-MAC® S IMAGER 8402 XS/8403 XS is not compatible with the C1200E General Processing Container/Tray.
-  **NOTE:** The Quick Connect Kit is not required as the instrument does not have any channels.

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6.11 Reprocessing limits

The product's service life and correct functioning are largely determined by mechanical stress and chemical influences within the scope of reprocessing and application.

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7 Service and repair

7.1 Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original parts.

7.2 Limitation of liability

KARL STORZ SE & Co. KG shall be liable for the failure or a deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the conditions that all assembly, operations, system expansions, readjustments, modifications, or repairs to the same have been performed by persons authorized by KARL STORZ, and that the instrument has been used in accordance with its operating instructions at all times.

7.3 Warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ. The medical device must always be sent to your local subsidiary (see “Subsidiaries” section), even during the warranty period.
Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

7.4 Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the device's reliability and extending its useful service life.
Information regarding maintenance services can be obtained from your local representative or from the manufacturer.

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7.5 Directive compliance

According to Medical Device Directive (MDD):

This medical device bears the CE mark in accordance with the Medical Device Directive (MDD) 93/42/EEC.

7.6 Standard compliance

(for 8402 XS/8403 XS)

According to IEC 60601-1, 60601-2-18,

UL 60601-1, CAN/CSA C22.2 No. 601.1-M90:

- Type of protection against electric shocks: N/A
- Degree of protection against electric shocks:
Applied part type BF

According to IEC 60601-1-2:

Please read the Electromagnetic Compatibility Information in the appendix.

Additional equipment connected to medical electrical equipment must be able to prove its compliance with the respective IEC or ISO standards (e.g., IEC 60950 for data processing equipment). Furthermore, all configurations shall comply with the standard requirements for medical systems (see IEC 60601-1-1 or Clause 16 of the 3rd Ed. of IEC 60601-1, respectively). Anybody connecting additional equipment to medical electrical equipment is a system configurator and is therefore responsible for the system's compliance with the standard requirements for systems. Please note that local laws take priority over the above-mentioned standard requirements. If in doubt, please consult your local specialist dealer or the technical service department.

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
7.7 Disposal

This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the device as electronic scrap. Please ask either KARL STORZ SE & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ SE & Co. KG is responsible for the proper disposal of this device.

7.8 Repair program

During the repair period, you will generally receive a unit on loan, which must then be returned to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to:

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Straße 34
78532 Tuttlingen, Germany

Outside of Germany, please contact your local KARL STORZ subsidiary or authorized dealer.

7.9 Important information

For hygiene reasons or to prevent infection, telescopes, instruments and devices must be cleaned, disinfected and, if necessary, sterilized before they are sent for repair.

We reserve the right to return contaminated instruments/devices to the sender.

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8 Technical description

C-MAC® S IMAGER 8402 XS

Classification:	IEC Class N/A
Power frequency:	DC
Voltage:	5V
Power consumption:	2W
Immersion protection:	IPX8
Illumination:	LED, white, 1 W
Camera technology:	CMOS
Resolution:	320 × 240 pixels
Weight:	225 g

C-MAC® S IMAGER 8403 XS

Classification:	IEC Class N/A
Power frequency:	DC
Voltage:	5V
Power consumption:	1 W
Immersion protection:	IPX8
Illumination:	LED, white, 0,2 W
Camera technology:	CMOS
Resolution:	300 × 400 pixels
Weight:	225 g

Operating/storage conditions for the C-MAC® S IMAGER

Operating conditions: 0 °C to +40 °C
30 % to 70 % rel. humidity

Storage conditions: -10 °C to +60 °C
5 % to 95 % rel. humidity

Barometric pressure: 700 – 1080 hPa

Average life cycle:

The end of the product's service life is largely determined by wear, reprocessing methods, the chemicals used and any damage resulting from use, so that a life cycle of up to 2 years can be reached.

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9 Appendix

9.1 Troubleshooting

WARNING: Always disconnect the device from the mains before carrying out any maintenance work!

Symptom:

Cloudy picture, stripes, streaks, or similar.

Possible causes:

The lens of the C-MAC® S video laryngoscope blade for single use is soiled.

The lens of the C-MAC® S IMAGER is soiled.

Remedy:

Use a new C-MAC® S video laryngoscope blade for single use.

Clean the optical surface of the image chip and LED illuminator on the C-MAC® S IMAGER with a cotton tip applicator soaked in 70 % isopropyl alcohol.

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9.2 Accessories

Art.no. 809125

MAGILL forceps, modified by BOEDEKER, length 25 cm, suitable for the endoscopic removal of foreign bodies, for use with video laryngoscopes of sizes 2 – 4

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9.3 Information on electromagnetic compatibility (EMC)

WARNING: Medical electrical devices are subject to special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this Appendix during installation and operation.

The C-MAC® S video laryngoscope blade for single use 051113-10/051114-10/051116-10 together with the C-MAC® S IMAGER 8402 XS/8403 XS complies with EN/IEC 60601-1-2 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The above-mentioned devices/combinations are Group 1 devices (acc. to CISPR 11).

Group 1 contains all the 'equipment and systems which use or generate RF energy exclusively for internal functions'.

NOTE: The tables and guidelines that are included in this Appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this device does cause electromagnetic interferences with other devices, the user is encouraged to try to correct the interferences by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the devices to different electrical circuits.

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WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: Devices must not be stacked on top of each other or operated close to each other, as this may cause EMC interference.

WARNING: The use of accessories other than those specified in this instruction manual may result in increased emissions or decreased interference immunity of the C-MAC® monitor 8403 ZX in combination with the C-MAC® S IMAGER 8403 XS. The accessories listed in this instruction manual have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories other than those specified here, it is the responsibility of the operator to ensure that they comply with IEC 60601-1-2.

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Guidance and manufacturer's declaration – electromagnetic emissions		
The C-MAC® S IMAGER 8402 XS/8403 XS in conjunction with the video laryngoscope blade for single use 051113-10/051114-10/051116-10 is intended for use in the environment specified below. The customer or user of the system should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Complies Group 1	The aforementioned devices and combinations of devices only use HF energy for their internal functioning. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Complies Class B	The aforementioned devices and combinations of devices are suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A. Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 2
Guidance and manufacturer's declaration – electromagnetic immunity

The C-MAC® S IMAGER 8402 XS/8403 XS in conjunction with the video laryngoscope blade for single use 051113-10/051114-10/051116-10 is intended for use in the electromagnetic environment specified below. The customer or user of the system should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment. If image distortion occurs it may be necessary to further filter the power quality by using a AC line filter.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T^* (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Complies <5 % U_T^* (>95 % dip in U_T) for 0.5 cycle Complies 40 % U_T (60 % dip in U_T) for 5 cycles Complies 70 % U_T (30 % dip in U_T) for 25 cycles Complies <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

* Note: U_T is the a.c. mains voltage prior to application of the test level.

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Table 4
Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting

The C-MAC® S IMAGER 8402 XS/8403 XS in conjunction with the video laryngoscope blade for single use 051113-10/051114-10/051116-10 is intended for use in the electromagnetic environment specified below. The customer or user should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the aforementioned devices or combinations of devices, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	$d = [3.5 / 3] \sqrt{P}$ 150 kHz to < 80 MHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = [3.5 / 3] \sqrt{P}$ 80 MHz to 800 MHz $d = [7 / 3] \sqrt{P}$ 800 MHz to 2,5 GHz
			where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

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Table 6
Recommended separation distances between portable and mobile RF communications equipment and the KARL STORZ C-MAC® S IMAGER 8402 XS/8403 XS in combination with the video laryngoscope blade for single use 051113-10/051114-10/051116-10

The C-MAC® IMAGER 8402 XS/8403 XS in conjunction with the video laryngoscope blade for single use 051113-10/051114-10/051116-10 is intended for use in an electromagnetic environment in which the HF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3.5}{3} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.1	0.1	0.2
0.1	0.4	0.4	0.7
1	1.2	1.2	2.3
10	3.7	3.7	7.4

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters [m] can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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KARL STORZ Endoscopy Canada Ltd.
1717 Millcreek Drive, Mississauga
Ontario, L5N 3R3, Canada
Phone: +1 905 816-4500, Fax: +1 905 858-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info.canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8526
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: info@kseaa.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopy Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@kseala.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantico, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorzcuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Lago Constanza No 326,
Col. Granada Del. Miguel Hidalgo
C.P. 11520 México D.F.
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: info@karlstorz.com.br

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA – Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Østensjøveien 15B
0661 Oslo, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
12739 Skärholmen, Sweden
14108 Kungens Kurva, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi OY
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)96824775
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestucio g. 59
08124 Vilnius, Lithuania
Phone: +370 68567000
E-Mail: LT-LV-info@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopie Nederland B. V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N. V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopy France S. A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing@karlstorz.fr

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S. A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Adria Eos d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Verovškova c. 60A
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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KARL STORZ ENDOSKOPE
East Mediterranean and Gulf S.A.L.
Block M, 3rd Floor
Beirut Souks, Weygand Street
2012 3301 Beirut, Lebanon
Phone: +961 1 999390, Fax: +961 1 999391
E-Mail: info@karlstorz-emg.com

KARL STORZ Endoscopy (South Africa) (Pty) Ltd.
P.O. 6061, Roggebaai 8012
Cape Town, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ Endoscopy Kasachstan
6, Saryarka str. BC „Arman“, office314
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai – United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, MPlaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 8 3823 8000, Fax: +84 8 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy Australia Pty. Ltd.
15 Orion Road Lane Cove NSW 2066
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd., Hong Kong
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 Hing Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd.
Beijing Branch,
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
Beijing, 100022, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shanghai Branch
36F, Tower 1 Grand Gateway
No. 1 Hongguiao Road, Xuhui District,
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Chengdu Branch
Room 805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
Chengdu, 6100414, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Shenyang Branch,
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
Shenyang, 110014, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd.
Guangzhou Branch
Room 02B, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
Guangzhou, 510620, People's Republic of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com.sg

KARL STORZ Endoscopy Korea Co. Ltd.
I-Park Bundang 1, 102dong, 8F
239 Jeongjail-Ro, Bundang-Gu
Seongnam, 13556, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: mingwang@karlstorz.tw

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: , Fax: +63 2 31745-00
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

KARL STORZ SE & Co. KG

Dr.-Karl-Storz-Straße 34
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.com
Web: www.karlstorz.com